




Maximum productivity. Trusted results.

Thermo Scientific TraceFinder Software


Maximum productivity. Trusted results.

Designed for fast and flexible compound screening and quantitation, Thermo Scientific™ TraceFinder™ software provides unique features and capabilities to support a wide range of applications and analytical laboratories.


Built with a customizable user interface, flexible method templates, comprehensive compound database, and access to extensive spectral fragmentation libraries, TraceFinder software allows operators of all experience levels to confidently drive laboratory productivity.




Thermo Scientific™ Orbitrap Tribid™ mass spectrometer



Thermo Scientific™ Orbitrap Exploris™ mass spectrometer series




Thermo Scientific™ TriPlus™ RSH autosampler with Thermo Scientific™ TSQ™ Triple Quadrupole mass spectrometer series



Thermo Scientific™ Q Exactive™ HF hybrid quadrupole-Orbitrap™ mass spectrometer series


Harmonize
your lab
productivity




Direct control for a wide range of systems

- Multi-channel LC
- Single Quad MS
- Triple Quad MS
- High-Resolution LC-MS
- Ion chromatography
- Gas chromatography
- High-Resolution GC-MS


Operating procedures and
training can be streamlined and
standardized, minimizing training
requirements and simplifying
lab operations.




Thermo Scientific™ Vanquish™ Core HPLC system
with Thermo Scientific™ TSQ Fortis™ Triple-Stage
Quadrupole mass spectrometer



Thermo Scientific™ ISQ™
EC/EM Single Quadrupole
mass spectrometer



Thermo Scientific™ Dionex™ ICS-6000 HPIC system with
Tablet Interface and Core with Thermo Scientific™
Orbitrap Exploris™ 120 mass spectrometer



Thermo Scientific™ Q Exactive™ GC Orbitrap™ GC-MS/MS system with
Thermo Scientific™ TRACE™ 1300 Series GC injectors

Improve data quality and increase confidence for any screening assay

Your targeted screening list may be long—but method setup doesn't have to be

Large lists of target compounds can be easily and accurately managed with the use of TraceFinder customizable compound databases. Target compounds can be added, removed, or edited within the database and subsequently included a new or existing processing method. Compound databases within TraceFinder software are designed to include all the information necessary to accurately define your target compounds and can be easily populated manually as well as through file import.

Formatted text files as well as instrument method export lists can be readily imported to populate a compound database, while non-formatted text files can easily be imported via the built-in mapping tool; eliminating the need for file modification or reformatting prior to import.

Database mapping/import tool

Mass List Mapping Tool

Template: **import** Refresh Preview

Separator Options: Comma Tab Semicolon Space Other

Data starts on row:

Preview

<Ignore Column>	<Ignore Column>	<Ignore Column>	<Ignore Column>	<Ignore Column>	<Ignore Column>
1_A_nFEN	T1: 219.152->55.373	219.152	55.373	219.152	1_A_nFEN
1_A_nFEN_13C6	T1: 225.152->84.208	225.152	84.208	225.152	None
1_BUP	T1: 468.310->396.210	468.31	396.21	468.31	None
1_BUP_D4	T1: 472.400->400.220	472.4	400.22	472.4	None
1_NBUP	T1: 414.274->211.040	414.274	211.04	414.274	None
1_NBUP_D3	T1: 417.243->83.262	417.243	83.262	417.243	None

Save Template and Return Cancel

Database view example

Compound Database - EFS_HRAM_Compound_Database


Compound Name	Peak Label	Peak Workflow	Associated Target Peak	Chemical Formula	MS Order	Precursor m/z	Product m/z	m/z	Adduct
640 Amphetamine	T1: 206.15394	TargetPeak	*	C13H19NO	ms1	0.000	206.15394	0.000	Hydrogen
641 Amphetamine	T1F3: 206.15394->160.12020	Fragment	*	C13H19NO	ms2	206.15394	160.12020	160.12020	Hydrogen
642 Amphetamine	T1F3: 206.15394->133.06482	Fragment	*	C13H19NO	ms2	206.15394	133.06482	133.06482	Hydrogen
643 Amphetamine	T1F3: 206.15394->105.07032	Fragment	*	C13H19NO	ms2	206.15394	105.07032	105.07032	Hydrogen
644 Amphetamine	T1F4: 206.15394->86.09697	Fragment	*	C13H19NO	ms2	206.15394	86.09697	86.09697	Hydrogen
645 Amphetamine	T1F3: 206.15394->74.09708	Fragment	*	C13H19NO	ms2	206.15394	74.09708	74.09708	Hydrogen
646 Amphetamine	T1F3: 206.15394->72.08145	Fragment	*	C13H19NO	ms2	206.15394	72.08145	72.08145	Hydrogen
647 Amphetamine	T1: 242.16115	TargetPeak	*	C10H19NO2	ms1	0.000	242.16115	0.000	Hydrogen
648 Amphetamine	T1F2: 242.16115->143.08295	Fragment	*	C10H19NO2	ms2	242.16115	143.08295	143.08295	Hydrogen
649 Amphetamine	T1F2: 242.16115->113.07111	Fragment	*	C10H19NO2	ms2	242.16115	113.07111	113.07111	Hydrogen
650 Amphetamine	T1: 338.21146	TargetPeak	*	C22H27NO2	ms1	0.000	338.21146	0.000	Hydrogen
651 Amphetamine	T1F1: 338.21146->193.10091	Fragment	*	C22H27NO2	ms2	338.21146	193.10091	193.10091	Hydrogen
652 Amphetamine	T1F3: 338.21146->176.0772	Fragment	*	C22H27NO2	ms2	338.21146	176.0772	176.0772	Hydrogen
653 Amphetamine	T1F3: 338.21146->115.05416	Fragment	*	C22H27NO2	ms2	338.21146	115.05416	115.05416	Hydrogen
654 Amphetamine	T1: 209.12845	TargetPeak	*	C11H19NO2	ms1	0.000	209.12845	0.000	Hydrogen
655 Amphetamine	T1F1: 209.12845->152.10728	Fragment	*	C11H19NO2	ms2	209.12845	152.10728	152.10728	Hydrogen
656 Amphetamine	T1F2: 209.12845->137.08382	Fragment	*	C11H19NO2	ms2	209.12845	137.08382	137.08382	Hydrogen
657 Amphetamine	T1F3: 209.12845->122.06023	Fragment	*	C11H19NO2	ms2	209.12845	122.06023	122.06023	Hydrogen
658 Amphetamine	T1: 214.03778	TargetPeak	*	C8H13NO2	ms1	0.000	214.03778	0.000	Hydrogen
659 Amphetamine	T1F1: 214.03778->168.03253	Fragment	*	C8H13NO2	ms2	214.03778	168.03253	168.03253	Hydrogen
660 Amphetamine	T1F2: 214.03778->100.99018	Fragment	*	C8H13NO2	ms2	214.03778	100.99018	100.99018	Hydrogen
661 Amphetamine	T1F3: 214.03778->68.04864	Fragment	*	C8H13NO2	ms2	214.03778	68.04864	68.04864	Hydrogen
662 Amphetamine	T1: 204.95771	TargetPeak	*	C8H13NO2	ms1	0.000	204.95771	0.000	Hydrogen-Loss
663 Amphetamine	T1F1: 204.95771->168.96855	Fragment	*	C8H13NO2	ms2	204.95771	168.96855	168.96855	Hydrogen-Loss
664 Amphetamine	T1F2: 204.95771->134.96127	Fragment	*	C8H13NO2	ms2	204.95771	134.96127	134.96127	Hydrogen-Loss
665 Amphetamine	T1F3: 204.95771->88.01546	Fragment	*	C8H13NO2	ms2	204.95771	88.01546	88.01546	Hydrogen-Loss

An example of the Compound Database, showing the range of information stored for each compound such as transitions, precursor and product ion targets, and adducts. Additional information can be applied to aid filtering and reporting, such as compound type and compound grouping.

Importing target lists from a range of sources is simplified through the use of TraceFinder's Mass List Mapping Tool.

- **Build**—Targeted screening methods can be built from multiple compound databases to easily include hundreds of compounds
- **Optimize**—TraceFinder target screening methods focus around four main search features, *m/z*, retention time, fragment ions, and isotopic pattern
- **Analyze**—Each search feature within the method can be easily viewed and modified to provide the right level of search requirements for any assay

Quickly and easily build targeted screening methods to monitor hundreds of compounds.



Method View - Targeted_Screening Example

Peak Filter Settings

Use RT Limits Search from minutes to minutes

Use Matrix Blank Amplifier

Chromatogram View Width minutes

Use Source CID Scans

Show all compounds

Unknown Screening

Include Unknown Screening

Target Screening Settings

Enabled	Database Name	
<input type="checkbox"/>	Clin_Tox_Endura_SRM	open
<input type="checkbox"/>	Clin_Tox_Quantiva_SRM	open
<input type="checkbox"/>	DefaultGC	open
<input type="checkbox"/>	DefaultLC	open
<input checked="" type="checkbox"/>	EFS_Database	open
<input checked="" type="checkbox"/>	EFS_HRAM_Compound_Database	open
<input type="checkbox"/>	GCMSMS Pesticide Analyzer 1001	open
<input type="checkbox"/>	Metabolite_Database	open
<input type="checkbox"/>	Toxicology_HRAM_Compound_Database_v1	open

Identification and Confirmation Settings

Peaks *m/z* Threshold Override

S/N Ratio Threshold

Mass Tolerance: ppm

Retention Time Identify Confirm Ignore if Not Defined

Window Override (sec)

Fragment Ions Identify Confirm Ignore if Not Defined

Min. # of Fragments

Intensity Threshold

Mass tolerance ppm

MS Order

Isotopic Pattern Identify Confirm Fit Threshold (%)

Allowed Mass Deviation (ppm)

Allowed Intensity Deviation (%)

Use Internal Mass Calibration

Library Search Identify Confirm Library Search Type:

Targeted screening using multiple compound databases is quick and easy, with full flexibility to control identification and additional confirmatory settings such as isotopic pattern and library searching.

Harness the full power of high-resolution accurate-mass spectrometry with mzVault library searching


- High-resolution accurate-mass (HRAM) mass spectrometry is a powerful tool for any screening assay, and when coupled to the highly curated library content of **Thermo Scientific™ mzVault™** library you can build unprecedented confidence in your screening results.
- Access the MS²-level content of **Thermo Scientific™ mzCloud™** mass spectral fragmentation library offline, through mzVault
- Easily search multiple mzVault libraries for any screening method
- Generate custom mzVault libraries, enabling new and expanded compound search capabilities to support novel research

Application - Library Selections

Libraries:

	Enabled	Library Name	Library Type
1	<input type="checkbox"/>	MAINLIB	NIST
2	<input type="checkbox"/>	EFS_HRAM_Spectra_Library.db	mzVault
3	<input type="checkbox"/>	EFS_Library.db	mzVault
4	<input type="checkbox"/>	LibraryManager_Toxicology_combined_v1.db	mzVault
5	<input type="checkbox"/>	mzVault_May2018_ver5.db	mzVault

The ability to search multiple libraries from NIST, through to offline mzCloud contents with mzVault or your own libraries adds confidence to any screening assay.




View fragmentation spectra and associated metadata, search, maintain or build your own proprietary libraries to utilize and share your knowledge.

Exploring the unknowns— in any workflow


The ability to identify features of interest in a single sample or across a large set of samples requires a wide-ranging set of tools. For unknown screening, TraceFinder has combined best in-class library search capabilities with a comprehensive set of data visualization tools to provide a reliable, straightforward workflow for unknown component detection.

Unknown component investigation is enhanced by a wide array of data visualization tools such as:

- Heat map tools
- Cross sample peak overlays
- Cross sample peak lists
- Isotope matching display
- Group averages



Through customizable options to visualize and interrogate large data sets, such as heat maps for group analysis, user defined retention time alignment, and simultaneous library searching, the process of identification and decision-making is simplified.




Data visualization tools are easily implemented to highlight features of interest and aid with putative identification. In addition to heat map displays, all identified components can be displayed in a single chromatographic view, sample groups can be averaged and compared to experimental controls, and fragmentation spectra can be searched against multiple mzVault spectral libraries.

Customizable quantitation for maximum results


Comprehensive quantitation tools designed to meet any workflow

TraceFinder software simplifies quantitation across LC, IC and GC-quadrupole and high-resolution mass spectrometers for large and small target lists through:

- Template-driven workflows
- Compound databases for method building
- In-depth data visualization capabilities
- Market specific terminology built into each workflow



Whether quantifying using data from triple quadrupole or high-resolution based mass spectrometers, data review is streamlined through flagging and data quality checks to allow rapid data interrogation and report generation.



Build your workflow to fit the application with a customizable user interface

TraceFinder software provides a wide range of features and capabilities, but not all features are required for all experiments. With an application-wide customizable user interface, the workspaces and data display can be configured to visualize as few or as many of the application features as necessary for the analyst and experiment type.


Utilizing a fit for purpose display greatly increases the productivity during system setup, training and SOP development, and assay data review.

Primarily focusing on targeted and untargeted screening workflows for veterinary drugs, but also screening for other toxicants and contaminants, TraceFinder software has allowed our lab to perform these workflows in a simple and efficient way; we have compound databases allowing us to perform fragment matching and confident compound identification without processing raw data files through complicated workflows across multiple software package.


—Laura Burns, Diagnostic Associate II, Iowa State University
Veterinary Diagnostic Laboratory

Full display

Focused display



The full view within the navigation tabs (left column) provide access to all functionality within the application, allowing full flexibility and control for method modification, data acquisition, and data visualization.



The customized or focused view provides the capability to display only the required navigation tabs (left column), providing a fit for purpose display for increased operational simplicity.

Combine quantitation and targeted screening with TraceFinder software

With TraceFinder software a targeted screening workflow can be built directly into a general quantitation method providing ultimate assay versatility. Target compounds are easily selected from any available TraceFinder compound database creating an easy to view list for all compounds to be screened and quantitated in the assay. Target compounds can be added, removed, or edited all within the processing method.

Screening assays built with quantitation allow for full and single point calibration curves, customizable data review tools and compound flagging.

The screenshot displays the TraceFinder software interface. On the left, there is a 'Compound Database' pane with a tree view of compounds. The main area shows a 'Peak List' table with columns for Peak Label, Peak Workflow, Associated Target Peak, Chemical Formula, MS Order, Precursor m/z, Product m/z, Adduct, and Polarity. Below the table, there is a 'Compound Details' pane for the selected compound, Acetophenone, showing its chemical formula, CAS No., and compound type.

Target compounds are easily selected from any available TraceFinder compound database, creating an easy-to-view list of all compounds to be screened and quantitated in the assay. Target compounds can be added, removed, or edited all within the processing method, and all necessary targeting information is contained within the database.



Navigating within TraceFinder software is intuitive. It is a straightforward process to create master methods and process data. It is very convenient to be able to associate information from a compound database into a master method. The data processing allows for viewing data in different ways all in one screen.

—Dwayne Schrunck, Laboratory Director, Iowa State University Veterinary Diagnostic Laboratory

The screenshot displays the TraceFinder software interface for data review. The top pane shows a 'Sample Results' table with columns for Compound, RT, Type, Sample Results, Status, Confirm, PK, Flags, Sample Type, Filename, Comments, Excluded, Area, Respon, Actual RT, % CV, Level, %Diff, %RSD, and Calculated Amt. Below the table, there are three plots: a 'Data Review' plot showing a peak at RT 2.68, a 'Calibration Curve' plot showing a linear relationship between area and concentration, and a 'Peak List' table for the selected compound, Acetophenone.


Combining targeted screening and quantitation workflows provides both compound flagging for rapid compound review and calibration curve generation in a single processing method and data review workflow.




Customizable reporting for crucial results

The ability to easily and accurately report assay results is as essential as the data acquisition itself and is often a laboratory bottleneck. In addition to standard reporting templates, customized report generation is an integrated component within TraceFinder software saving time and resources, minimizing the potential for error, and maximizing throughput.

- Integrated custom report generation saves time and resources, minimizing the potential of error, and maximizing throughput
- A broad selection of pre-configured reports are available and ready for immediate use or can be used to provide a starting template for further report customization



With intuitive and predefined options for inserting chromatograms and building data tables, custom reports can be easily built to meet a wide range of reporting requirements. A broad selection of pre-built reports is available and ready for immediate use or can be used to provide a starting template for further report customization.




Report generation provides an additional tool for assay evaluation that can be built into the overall data review process. Using additional reporting features such as custom logic statements, data filters, and custom calculations, the custom reporting capability can provide a high-level overview of assay results and to identify samples that require additional interrogation all within the data review workspace.

Maximum productivity. Trusted results.

Maximize your connectivity with TraceFinder software and B-Link


TraceFinder software not only has the ability to control multiple inlets (IC, LC, GC and multi-channel LC) and types of instrumentation (single quadrupole, triple quadrupole and high-resolution MS), it also allows bi-directional flow of methods and results from LIS/LIMS (Laboratory Information System/Laboratory Information Management System) to further enhance your productivity.

- Quantitative methods and screening
- Transfer of results and quality control values
- Minimizing user settings and interactions
- Simple to install, maintain and to use



B-link (BYG INFORMATIQUE, L'UNION France), a universal instrument connector, enables TraceFinder software to connect any LIS, LIMS or middleware, allowing the transfer of requests, information, and results. These requests can be in real time or in batch mode, without any user-action required within B-link once an initial simple system integration has been performed.

Connectivity between B-Link Connector and TraceFinder software



Challenge

Communication between the MS platform and LIS/LIMS

- Downloading test requests
- Uploading test reports

B-Link LIS/LIMS Connector

Software providing bi-directional communication

- Between LIS/LIMS and B-Link
- Between B-Link and TraceFinder software

Don't wait,
try TraceFinder
software
today.

Whether you are a new user, or looking to upgrade, try the latest version of TraceFinder software for free. With your **free 60-day trial**, access TraceFinder software's full functionality to see what knowledge and insights you can discover, and how the latest productivity improvements can help improve laboratory productivity and drive confidence.

Download it today: thermofisher.com/TraceFinder





Almanac web-based monitoring and management

Stay connected to your science. See how the Thermo Scientific™ Almanac™ application can help you get the most out of your instruments.

thermofisher.com/almanac



Technical and online support: peak performance for your instruments

Helping you keep your instruments running at peak performance is our goal. Whether you're looking for an instrument manual or spare parts, want to submit a repair request, or check on the status of your warranty or service contract, we have every support option you're looking for.

thermofisher.com/technicalresources



Protecting your investments: unparalleled laboratory services

Unity™ Lab Services provides a single source for integrated lab service, support, and supply management. Our customized service offerings and world-class service experts have the flexibility and experience to address your laboratory's needs. We provide a complete portfolio of services and support solutions designed to help you improve productivity, reduce total cost of ownership, and ensure performance throughout your laboratory.

unitylabservices.com

Thermo Scientific™ Compound Discoverer™ software

thermofisher.com/CompoundDiscoverer

Thermo Scientific™ Mass Frontier™ software—Library curation capabilities

thermofisher.com/MassFrontier

Thermo Scientific™ mzVault Offline Mass Spectral database

thermofisher.com/mzVault

Thermo Scientific™ Chromeleon™ CDS & compliance-ready screening/quantification

thermofisher.com/Chromeleon

Find out more at thermofisher.com/TraceFinder