



# GRABBER D885

Enhance your PAL RTC  
Grabb anything which is not magnetic

## Features

- Special tool designed for all PAL RTC
- Compatible with GC and LC technique
- Transport of non-magnetic objects
- Provided with standard tasks
- Customized method available

## Examples of application:

- Automatic filtering of solutions using Thomson vials
- Transport of test tubes
- Transport of Spectroscopy cuvettes
- Transport of standard caps vials
- and much more...



You need a special tool ?

Grabber D885 can be configured to all your needs

Load a sample, filter it and inject it in a single method

Delivered with the following Pal Sample Control activities:

- Filter
- Grab and Drop
- Non-magnetic transport
- Easy customization of the tool to most shape
- Up to 40mm wide objects: **Challenge us!**



palsystem@brechbuehler.ch www.brechbuehler.ch

**BRECHBÜHLER**  
AG

a Value Added Reseller of

**PAL** SYSTEM  
Ingenious sample handling

SCIENTIFIC ANALYTICAL SOLUTIONS